

HCP

	HCP exkl KM			HCP exkl KM			HCP exkl KM			HCP exkl KM			HCP exkl KM							
	2025 VT 70% diff 185			2025 VT 70% diff 185			2025 HT 70% diff 185			2025 HT 70% diff 185			2026 VT 70% diff 185			2026 VT 70% diff 185				
	3 s	1 s	Jan/Feb	3 s	1 s	Mar/Apr	3 s	1 s	Sep/okt	3 s	1 s	Nov/Dec	3 s	1 s	Jan/Feb	3 s	1 s	mar/apr		
Agneta Bergqvist																				
Anita Långström	481	25	18	477	26	18	461	31	22	442	38	27	450	35	25	426	43	30		
Ann Rosenquist	429	42	30	428	42	30	441	38	27	382	58	41	396	53	37	414	47	33		
Ann Sundqvist/Sundelin	475	27	19	496	20	14	487	23	16	432	41	29	450	35	25	492	21	15		
Anneli Vestman	323	78	45	299	85	45	312	81	45	305	83	45	300	85	45	0	0	0		
Ann-Marie Johander	449	35	25	468	29	20	459	32	22	459	32	22	450	35	25	467	29	20		
Ann-Sofie Asplund	549	2	1	574	-6	1	631	-25	1	518	12	8	519	12	8	536	6	4		
Birgitta Pettersson	413	47	33	326	76	45	337	73	45	347	69	45	357	66	45	373	61	43		
Catrine Cederberg	493	21	15	472	28	19	463	31	21	421	45	32	435	40	28	446	36	25		
Eva Palm Dahlström	501	18	13	533	7	5	0	0	23	500	18	13	510	15	11	494	20	14		
Gerd Lundqvist	472	28	19	447	36	25	447	36	25	448	36	25	453	34	24	483	24	17		
Gunilla Törnfeldt	50%	0	0	23	50%	410	49	33	420	45	32	421	45	32	417	46	32	421	45	32
Inger Enström		462	31	22		487	23	16	479	25	18	471	28	20	474	27	19	487	23	16
Kristina Zettergren	Ny	240	105	45		257	99	45	279	92	45	269	95	45	270	95	45	316	80	45
Lena E Johnsson		360	65	45		349	69	45	336	73	45	328	76	45	333	74	45	339	72	45
Lena Gullstrand		286	90	45		325	77	45	297	86	45	334	74	45	339	72	45	314	80	45
Lena Thurn		492	21	15		512	14	10	517	13	9	488	22	15	486	23	16	495	20	14
Mari Nilsson		344	71	45		403	51	36	394	54	38	332	74	45	348	70	45	378	59	41
Mari Svensson		398	52	37		409	49	34	395	53	37	433	41	29	426	43	30	410	48	34
Monika Eriksson		0	0	0		399	52	36	0	0	23	431	41	29	423	44	31	420	45	32
Ulla Sundqvist/Säkvist	50%	0	0	23	50%	483	24	17	477	26	18	455	33	23	459	32	22	434	40	28
Åsa Kolmodin													287	62	45	327	76	45		

Om färre matcher än hälften endast 50% av maxhcp (23).

0,7% på diff 185 blir det nya hcp, max 45 hcp